



WWW.SRITRANGGLOVES.COM



Company Profile **2019**  
Update APR 2019

The logo consists of the letters 'S', 'T', 'G', and 'T' in a bold, purple, sans-serif font. The 'S' is positioned at the top left, the 'T' is to its right and slightly lower. Below the 'S' is a 'G', and to its right is another 'T'. The letters are arranged in a staggered, overlapping manner. The background features a faint, light purple hexagonal pattern.

STGT

Specialist & Teamwork in Giving  
「Touch Of Life™」



「Touch Of Life™」

In life, every touch has a meaning.  
From the first day that we are born into this world,

---

A pair of gloves takes us to our mother's embrace.  
As you move through stages of life,  
gloves help to protect you. Across industries, from medical to food and manufacturing  
of many great products, gloves have become an essential part of our life.



SINCE 1989,

**Sri Trang Gloves (Thailand) Public Company Limited**

takes pride in our attention to detail at every production step;  
from the careful selection of raw materials to the constant research and  
innovation. We take pride in producing gloves that deliver 「Touch Of Life™」

Not just a glove manufacturer,  
we are also part of a

**“Fully integrated  
rubber company”**





OUR VISION IS TO DELIVER  
「Touch Of Life™」  
ON A GLOBAL SCALE

Since  
1989

4000+  
Employees

123  
Production  
Lines

21  
Billion  
Pieces Per  
Annual

365  
Million USD  
Turnover  
2018

Export to  
**120+**  
Countries

## STGT LOCATION



OUR FOCUS



# WE RESEARCH INTO THE SCIENTIFIC HEART OF RAW MATERIALS

We believe that great products all start with raw materials.

At **SRI TRANG GLOVES (THAILAND)**

we heavily invest in our R&D, from choosing the right type of rubber trees, the research and development of plantation techniques, which the aim of producing the best quality latex with low protein content to supply to our Gloves factories. Not just latex, we also focus on Nitrile, using rigorous sourcing and testing to ensure that all Nitrile butadiene rubber from our suppliers is up to STGT standards.



Plantation Technique

We choose the best clone and develop plantation techniques to improve yield and quality.



Protein Control

We are expert in controlling every aspect of latex specification through our expert research and development to ensure that we create the best formula for each type of gloves.



Rigorous Sourcing

We use our expertise and detailed knowledge in sourcing and testing to ensure that all Nitrile Butadiene raw material from our Suppliers meets STGT stringent standards.





## CUTTING EDGE TECHNOLOGY AND NON-STOP INNOVATION

With years of experience, **SRI TRANG GLOVES (THAILAND)**

engineering team is the expert in designing of state-of-the-art glove production process. As a result, we are able to come up with highly flexible and efficient production process to ensure high quality hand protection to satisfy all our customers on a global scale.



Machine master

We are the expert in producing highly flexible and efficient production process to ensure high quality hand protection.



Automation

We take pride in our development of the automation process from production to auto packing with the aim of improving our gloves, reducing defect and delivering the "Touch Of Life" products on a global scale.



Non-stop melioarate technology

We never stop to improve the process and are committed to continuous product improvement to improve the quality of hand protection and reduce costs in order for our "Touch Of Life" products to reach all users around the world.



## OUR PRIORITY IS TO DELIVER THE 「Touch Of Life™」 and sustainable world

We are committed to conducting business with transparency and fairness to produce quality and environmentally friendly products that deserve to be call the 「Touch Of Life™」 and we care about all stakeholders in the supply chain, from upstream to downstream.  
Not just for the people but for every life in this beautiful world.



### Sustainability for farmers and fair trade

At SRI TRANG GLOVES (THAILAND), we have high standards and expect our suppliers to comply with our on principles on human rights, environmental impact and business integrity, all part of our 「Touch Of Life™」



### The environmentally friendly technologies

Apart from production process, SRI TRANG GLOVES (THAILAND) also invests heavily in biomass energy and waste management and green energy to ensure positive environment within our facilities as well as communities around our facilities.



### Human rights

We see our employees as part of the SRI TRANG GLOVES (THAILAND) family. We are committed to eliminating child labor and other abuses where they are found, and to achieving safe and healthy working conditions as well as following the 「Touch Of Life™」 principles. In addition to implementing our human rights commitment, we also joined the Business Social Compliance Initiative (BSCI).



### Animal Welfare

At STGT, we appreciate and love every life that makes up this wonderful world. We strongly commit to delivering the 「Touch Of Life™」 to many of unfortunate animals around Thailand. This includes supporting charities on stray animals, wild life and we are proud to be part of a better world.



### Sharing Touch of Love

At STGT, apart from creating 「Touch Of Life™」 products to serve customers globally, we also take big part in sharing to those unfortunate ones.



**SRI TRANG**  
GLOVES THAILAND

## DESIGNING A HAND PROTECTION FUTURE

We offer a wide range of high-quality natural latex and nitrile examination gloves

as well as general-purpose disposable gloves, both powdered and powder-free, catering to a wide variety of applications, including healthcare and medical, food processing, food handling, manufacturing and operations where hygiene standards come first.



# LATEX POWDERED

Latex Powdered has been part of our life for more than 100 years.

It was one of the life-changing innovations that still deliver the touch of life to people around the world today. Today latex gloves have been improved compared to when they were first produced. STGT takes pride in our knowledge and R&D in latex raw materials as well as experience in gloves production. We deliver the most comfortable gloves that remain effective in order to carry on this legacy and offer the 「Touch Of Life™」 on a global scale.

## Our Standard Latex Powdered Gloves



### PREMIUM LATEX POWDERED

Through years of research, we select the most suitable Latex formula, which combines a perfect thickness while retaining high flexibility and comfort.

**Properties** PREMIUM COMFORT POWDERED  
 Average M size: **5.8G** | Weight **240mm** | FINGER TEXTURED

---

**Used**

---

**Standard** ASTM D3578 ISO 1197 EN 455 MDD CE PPE

CE 2777 | PPE



### BALANCE LATEX POWDERED

This flagship of the powdered gloves is suitable for all kind of industries, including full medical use as examination gloves in hospitals

**Properties** BALANCE COMFORT POWDERED  
 Average M size: **5.0G** | Weight **240mm** | FINGER TEXTURED

---

**Used**

---

**Standard** ASTM D3578 ISO 1193 EN 455 MDD CE PPE

CE 2777 | PPE



### COMFORT LATEX POWDERED

Our Comfort Latex PD is designed to mimic no less marvel than the human skin.

**Properties** FINGER TEXTURED HIGH COMFORT POWDERED  
HIGH ELASTICITY HIGH SENSITIVITY | Average M size **4.8G** | Weight **220mm** | Length **230mm**

---

**Used**

---

**Standard** ASTM D3578

For more of our 「Touch Of Life™」 Products please contact [glove@sritranggroup.com](mailto:glove@sritranggroup.com)

# LATEX POWDER – FREE

Latex Powder-Free are designed to be alternative part of life

changing products. They offer more comfortable Touch for user without leaving any residue. With improved technology, Latex Powder-free gloves offer smooth easy donning for inner surface finish.

## Our Standard Latex Powder-Free Gloves



### LOW FRICTION LATEX POWDER-FREE

If the main goal is speed and maximum protection, the Low Friction Latex PF solution is the name of the game.

**Properties**

EASY DONNING	LOW FRICTION		
Average M size	6.2G	240mm	FINGER TEXTURED

**Used**

**Standard**

ASTM D3578 | ISO 1193 | EN 455 | MDD | CE | PPE | UH | CE 2777 | PPE



### GRIP PLUS LATEX POWDER-FREE

An all-rounder, with a stylish white appearance and optimal grip textured surface, suitable for medical, protective and food contact applications.

**Properties**

MAXIMUM GRIP			
Average M size	5.8G	240mm	FULL TEXTURE

**Used**

**Standard**

ASTM D3578 | ISO 1193 | EN 455 | MDD | CE | PPE | UH | CE 2777 | PPE



### COMFORT LATEX POWDER-FREE

The unique Latex powder free glove, a masterpiece of know-how and production process control that fulfills the requirements of a medical examination and a category III protective glove, including safe food handling.

**Properties**

EASY DONNING	HIGH COMFORT		
Average M size	5.0G	240mm	FINGER TEXTURED

**Used**

**Standard**

ASTM D3578 | ISO 1193 | EN 455 | MDD | CE | PPE | UH | CE 2777 | PPE

For more of our 'Touch Of Life' Products please contact [glove@sritranggroup.com](mailto:glove@sritranggroup.com)

## NITRILE POWDER-FREE

Apart from Latex gloves,

SRI TRANG GLOVES (THAILAND) has been producing Nitrile gloves,

using rigorous sourcing and testing to ensure that every Nitrile from our suppliers is up to STGT standards.

Our Nitrile examination gloves offer added value to users in term of “comfort & elastic feeling” as well enhanced protection against chemicals, and have been tested for resistance to selective chemicals and chemotherapy drugs (ASTM F739)\*

## Our Standard Nitrile Powder-free Gloves



### FORCE NITRILE POWDER-FREE

The tank of disposable nitrile gloves, a powerful protection against all kinds of harm.

Properties	Average M size	Weight	Length	FINGER TEXTURED
	6.2G	240mm		
Used	EASY DONNING	HIGH VISIBILITY	SUPERIOR FORCE AT BREAK	
Standard	ASTM D6319	ISO 11193	EN 455	MDD CE PPE
	CE 2777   PPE			



### BALANCE NITRILE POWDER-FREE

This legendary product is famously used in laboratories worldwide. Its balance of protection and tactile sensitivity convinced users at all different kind of applications.

Properties	Average M size	Weight	Length	FINGER TEXTURED
	4.8G	240mm		
Used	EASY DONNING	HIGH VISIBILITY		
Standard	ASTM D6319	ISO 11193	EN 455	MDD CE PPE
	CE 2777   PPE			



### COMFORT NITRILE POWDER-FREE

Our Comfort Nitrile PF are ideal for those who love reliable protection with great skin alike comfort and reduced fatigue in your hand.

Properties	Average M size	Weight	Length	FINGER TEXTURED
	3.5G	240mm		
Used	EASY DONNING	HIGH VISIBILITY	HIGH COMFORT	
Standard	ASTM D6319	EN 455	MDD CE PPE	
	CE 2777   PPE			



### COMFORT PLUS NITRILE POWDER-FREE

This Comfort+ Nitrile PF is designed to let you forget you are wearing gloves.

Properties	Average M size	Weight	Length	FINGER TEXTURED
	3.0G	240mm		
Used	PREMIUM COMFORT	EASY DONNING	HIGH VISIBILITY	
Standard	ASTM D6319	EN 455	MDD CE PPE	
	CE 2777   PPE			

For more of our 'Touch Of Life' Products please contact [glove@sritrangroup.com](mailto:glove@sritrangroup.com)

## STERILE EXAMINATION GLOVES

**STERILE R**

$\gamma$  min 25 kGy  
Gamma

Hygiene and cleanliness of our sterile gloves are our first priority.

At STGT, we choose Gamma ray radiation process that has high penetration and has no residue left over in order to ensure that our sterilized gloves are clean from germs and chemicals. All our sterile gloves are processed and monitored by a global leading sterilization and contamination specialist with more than 40 years of experience in providing safe and effective gamma radiation sterilization.



### Packing of sterile gloves.

POLYOLEFIN film is used for packing each pair of our sterile gloves to ensure the hygiene and safety of each glove. Polyolefin is recognized to be the packaging material that can withstand the sterilization from ETO, Electron beam to high penetration Gamma ray.

# TABLE OF PRODUCTS

## SPECIFICATION

## USAGE

## STANDARD

Product	Material	Colour	Inner Treatment	Outer Surface	Length	Weight (M SIZE)	Thickness at palm (Single wall)		AQL Up to	Proteins content	Personal Protection	Hospital/ Dental	Laboratory / Pharma	Food / Catering	Cleaning	Beauty/ Body Art	MDD (medical device directive 90/269/EEC)	PPE (Personal protective equipment 89/656/EEC)	Food contact (EC No. 1831/2003)	Food contact (ISO 22000 Food safety management)	EN 455	EN 420	EN 374-2	EN 374-3	ASTM F1671	ASTM D3578	THAILAND GMP	ABNT NBR 13392		
							Specification	Actual Production Value																						
Premium Latex Powdered	Natural Latex	Natural Latex White	Powdered	Finger Textured	Median 240mm	5.8 ± 0.2g	0.08 mm 3.1 mil	Median 0.11 mm 4.3 mil	1.5	</= 200 ug/g	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a	
Balance Latex Powdered	Natural Latex	Natural Latex White	Powdered	Finger Textured	Median 240mm	5.0 ± 0.2g	0.08 mm 3.1 mil	Median 0.10 mm 3.9 mil	1.5	</= 200 ug/g	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Comfort Latex Powdered	Natural Latex	Natural Latex White	Powdered	Finger Textured	- Minimum 220 mm For XS, S - Minimum 220 mm For M, L, XL	4.8 ± 0.2g	0.08 mm 3.1 mil	Median 0.09 mm 3.5 mil	1.5	</= 200 ug/g	☑	☑	☑	☑	☑	Class 1	n/a	n/a	☑	n/a	n/a	n/a	n/a	☑	*	☑	☑	☑	n/a	
Low Friction Latex Powder-Free	Natural Latex	Natural Latex White	Polymer Coated + 2x Chlorination	Finger Textured	Median 240mm	6.2 ± 0.2g	0.08 mm 3.1 mil	Median 0.12 mm 4.7 mil	1.5	</= 50 ug/g	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Grip Plus Latex Powder-Free	Natural Latex	Natural Latex White	Polymer Coated + Chlorination	Full Textured	Median 240mm	5.8 ± 0.2g	0.08 mm 3.1 mil	Median 0.11 mm 4.3 mil	1.5	</= 50 ug/g	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Comfort Latex Powder-Free	Natural Latex	Natural Latex White	Chlorination	Finger Textured	Median 240mm	5.0 ± 0.2g	0.08 mm 3.1 mil	Median 0.10 mm 3.9 mil	1.5	</= 50 ug/g	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Force Nitrile Powder-Free	Nitrile	Blue	2 Sided Chlorination	Full Textured	Median 240mm	6.0 ± 0.2g	0.05 mm 2.0 mil	Median 0.14 mm 5.6 mil	1.5	Not Present	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Balance Nitrile Powder-Free	Nitrile	Blue	Chlorination	Finger Textured	Median 240mm	4.8 ± 0.2g	0.05 mm 2.0 mil	Median 0.09 mm 3.5 mil	1.5	Not Present	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Comfort Nitrile Powder-Free	Nitrile	White, Blue, Pink Lady, Ocean blue, Violet blue, Black	Chlorination	Finger Textured	Median 240mm	3.5 ± 0.2g	0.05 mm 2.0 mil	Median 0.08 mm 3.1 mil	1.5	Not Present	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑	n/a
Comfort Plus Nitrile Powder-Free	Nitrile	Dark blue, White	Chlorination	Finger Textured	Median 240mm	3.0 ± 0.2g	0.05 mm 2.0 mil	Median 0.07 mm 2.8 mil	1.5	Not Present	☑	☑	☑	☑	☑	Class 1	Category III Complex Design	☑	☑	n/a	☑	n/a	☑	☑	☑	☑	☑	☑	☑	n/a

\* test base on in-house method

BECAUSE WE ARE  
**A PART OF YOU**



「Touch Of Life™」



## STANDARDS

Different countries have different standards,  
in order to deliver our 「Touch Of Life™」 products on a global scale,  
we take a close look to make sure our products are compliant with all standards.



**EN  
455**



As the goal of delivering 「Touch Of Life™」 products on a global scale. We produce disposable gloves that are compliant to worldwide standards however we also produce gloves that are surpass these as well.

## Factory Management Certificates



ISO 9001:2015  
Quality Management System



DIN EN ISO 13485:2016  
EU Medical devices Quality Management System



ISO 13485:2016  
Medical devices Quality Management System



ISO 22000:2005  
Food Safety Management System Cover GMP+HACCP



ISO 17025:2005  
Laboratory Quality Standard for Medical device testing



GMP  
Good Manufacturing practices



PPE (EU) 2016/425  
Module D  
EC quality monitoring conformity to type based on quality assurance of the production process

## Product Certificates (Medical)



uon. 1056-2556  
TISI  
Thai Industrial Standards



MDD 93/42/EEC Class I  
European Standard for Medical device directive Class I



MDD 93/42/EEC Class I Sterile  
European Standard for Medical device directive Class I in Sterile conditions

## Product Certificates (Industrial)



PPE (EU) 2016/425 Category I  
European Standard for Personal Protective Equipment Category I



PPE (EU) 2016/425 Category III  
EU Type-Examination certificate shown to satisfy the applicable essential health and safety requirements

## Sustainability Certificates



ISO 14001:2015  
Environmental Management System



Business Social Compliance Initiative  
Supply Chain management system international labour standards protecting workers' rights (STGT receive Rank A)



Green Industry  
Sustainability development for Manufacturer network from Thai Ministry of Industry (STGT receive LEVEL 3)



ECOVADIS  
Business Sustainable Rating (STA group received Silver rating for Natural Rubber)



Environmental Good Governance  
Sustainable Business Certificate From Thai Ministry of Industry

# CERTIFICATES

With years of hard work and our expertise commitment

to creating a complete chain of sustainable gloves business, Sri Trang Gloves (Thailand) is proud to have received numerous certificates and awards from many differences fields.

# AWARDS

## SET

CSR Recognition  
Award Rising Star 2014



## SET

Sustainability  
Award 2015



## SET

Thailand Sustainability  
Investment award  
2018, 2017, 2016, 2015



## SET

Rising Star  
Sustainability  
Award 2018



## Prime Minister

Best Exporter award 2018

## Prime Minister

Industry award 2018



## CSR-DIW

Continuous  
Award 2018



## CSR-DIW

Award 2016



## THAI FDA

Quality award  
2018, 2017, 2016,  
2014, 2011, 2010, 2009



# CONTACT US

## Factory



THAILAND FACTORY

### **SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED (HAT YAI BRANCH)**

No. 110, Kanjanavanit Road, Pahtong, Hat Yai, Songkhla, 90230 Thailand  
Tel. 66-7447-1471, 66-7429-1648-9, 66-7429-1471-5  
Fax. 66-7429-1650

### **SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED (SURAT THANI BRANCH)**

No. 189, Moo 7, Phlai Wat, Kanchanadit, Surat Thani, 84160 Thailand  
Tel. 66-7727-7888  
Fax. 66-7727-7886, 66-7727-7889

### **SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED (TRANG BRANCH)**

No. 85, Moo 6, Khuan Thani, Kantang, Trang, 92110 Thailand  
Tel. 66-7520-1199  
Fax. 66-7520-1150



## Office



THAILAND OFFICE



USA OFFICE



CHINA OFFICE



VIETNAM OFFICE

### **SRI TRANG GLOVES (THAILAND) PUBLIC COMPANY LIMITED (BANGKOK BRANCH)**

17<sup>th</sup> Floor, Park Venture Ecoplex Unit 1701, 1707-1712  
57 Wireless Road, Lumpini, Pathumwan, Bangkok, 10330, Thailand  
Tel. (66) 0-2207-4500 Fax. (66) 0-2108-2241-44

### **Sri Trang USA, Inc.**

5401 West Kennedy Boulevard, Suite 760 Tampa, FL 33609-2447  
Tel. +1 (813) 606-4301 Fax. +1 (813) 606-4314

### **SHIDONG SHANGHAI MEDICAL EQUIPMENT CO., LTD.**

5F East, Block D&F, W square No.1686 Wuzhong Road,  
Minhang District, Shanghai, 201103  
Tel. 021-6413-7860-603 Fax. 021-6413-7315

### **SRI TRANG INDOCHINA (VIETNAM) COMPANY LIMITED (HO CHI MINH BRANCH)**

Room No. 7.01A, 7th Floor, Vietnam Business Center Building.  
57-59 Ho tung Mau Street, Ben Nghe Ward, District 1, HCMC, Vietnam  
Tel. +84 8 3821 6869 Fax. +84 8 3821 6877

To Discover our "Touch Of Life" products line Please visit [www.srtranggloves.com](http://www.srtranggloves.com) or contact [Glove@srtranggroup.com](mailto:Glove@srtranggroup.com)