Quick reference guide for acute care facilities

Cardinal Health™ Exam Gloves



Clinical staff safety is too critical to compromise

Healthcare workers within acute care facilities perform many different tasks, so they need a strategic portfolio of gloves to meet the needs of specific jobs to be performed, working environment, duration of use and potential hazards.¹ Cardinal Health recommends a product mix from two different categories: high use nitrile and specialty.

High use nitrile portfolio

Non-sterile exam gloves designed to serve a broad range of everyday clinical needs. Cardinal Health recommends selecting **two different products** within this category: one FLEXAL™ glove for general use and one ESTEEM™ glove for extra protection when needed.

Fingertip thickness (+/-0.2mil)	Product family	Size	Cat. no.	Features	Color	Texture/ cuff length	Chemo drug tested (ASTM D6978)	Chemical tested (ASTM F739)	Packaging
General use	•								
3.7 mil	FLEXAL™ Comfort	XS	88FS01XS		Light blue	1.5%	USP <800>	✓	XS – L: 250 ea/bx; 10 bx/cs
		S	88FS02S	Dispensing film >650% elongation Low dermatitis potential					
		М	88FS03M						
		L	88FS04L						XL - XXL:
		XL	88FS05XL						230 ea/bx;
		XXL	88FS06XXL						10 bx/cs
3.5 mil	FLEXAL™ Touch with film	XS	88RTF01XS	Dispensing film	Blue	15.5	USP <800>		XS – L:
		S	88RTF02S						250 ea/bx;
		М	88RTF03M						10 bx/cs
		L	88RTF04L						XL – XXL:
		XL	88RTF05XL						230 ea/bx;
		XXL	88RTF06XXL						10 bx/cs
3.5 mil	FLEXAL™ Touch	XS	88RT01XS		Blue	,5%	USP <800>		XS – L:
		S	88RT02S	1					250 ea/bx;
		М	88RT03M						10 bx/cs
		L	88RT04L						XL - XXL:
		XL	88RT05XL						230 ea/bx;
		XXL	88RT06XXL						10 bx/cs
Extra prote	ction				,				
4.5 mil	ESTEEM™ Comfort	XS	88ES01XS	Dispensing film >650% elongation Low dermatitis potential	Light blue	15g	USP 50.6 <800> min	✓	XS – L:
		S	88ES02S						200 ea/bx;
		М	88ES03M						10 bx/cs
		L	88ES04L						XL-XXL:
		XL	88ES05XL						180 ea/bx;
		XXL	88ES06XXL						10 bx/cs
4.5 mil	ESTEEM™ Stretch	XS	8854NXSB		Teal blue	, 28 , 28	USP <800>	✓	XS – L:
		S	8855NSB						150 ea/bx;
		М	8856NMB						10 bx/cs
		L	8857NLB						XL – XXL:
		XL	8858NXLB						130 ea/bx;
		XXL	8859NXXB						10 bx/cs

USP <800> compliant

Thiotepa minimum breakthrough detection time in minutes



Carmustine minimum Smooth breakthrough detection

















Specialty portfolio

Non-sterile and sterile exam gloves designed to meet protection needs for specific clinical applications. The most common clinical applications are chemotherapy and chemical protection, sterile procedures and decontamination.

Fingertip thickness (+/-0.2mil)	Clinical application	Product family	Size	Cat. no.	Color	Texture/ cuff length	Chemo drug tested (ASTM D6978)	Chemical tested (ASTM F739)	Packaging
Chemother	apy and chemical prote	ection							
7.1 mil			S	8896EC		HIS*	USP 46.9 54.8 min		S – XL: 100 ea/bx; 10 bx/cs
	Chemotherapy handling and administration	ESTEEM™ Chemotherapy	М	8897EC	Dark blue				
			L	8898EC	— Dark blue				
			XL	8899EC					
	Chemical protection with an extended cuff length	Extended Cuff	XS	N8850XPB		, ₁₂	USP <800>	✓	XS – XL: 100 ea/bx; 10 bx/cs
6.9 mil			S	N8851XPB					
			М	N8852XPB	Blue				
			L	N8853XPB					
			XL	N8854XPB					
4.7 mil	Chemical protection with a standard cuff length	ESTEEM™	XS	8895NB		.556	USP 32.2 <800> min	✓	XS – L: 150 ea/bx;
			S	8896NB					
			М	8897NB	Blue				10 bx/cs
			L	8898NB	blue				XL – XXL: 130 ea/bx; 10 bx/cs
			XL	8899NB					
			XXL	8899NXXB					
Sterile									
	Sterile procedures/ compounding	Sterile (Pairs)	S	N8830		, in the second	USP 60.0 30.3 min min	✓	S – XL: 40 pr/bx; 5 bx/cs
6.5 mil			М	N8831	Blue				
			L	N8832					
			XL	N8833					
6.5 mil	Sterile procedures/ labor and delivery exams	Sterile (Singles)	S	N8820		الله الله الله الله الله الله الله الله	USP 60.0 30.3 min	✓	S – XL: 50 ea/bx; 5 bx/cs
			М	N8821	Blue				
			L	N8822					
			XL	N8823	_				
Decontami	nation		7.2	1.10020				<u> </u>	ļ
11.0 mil	Nitrile glove for decontamination, tested for puncture resistance ³	Decontamination	S	88NDS		, <u>i</u>	USP 67.8 <800> min	✓	S-XL: 50 ea/bx; 10 bx/cs
			M	88NDM					
			L	88NDL	Blue				
			XL	88NDXL	_				
			S	88LDS		<u></u>			
18.0 mil	Latex glove for decontamination	Latex Decontamination⁴	M	88LDM	_		√		S – XL: 50 ea/bx; 10 bx/cs
			LG	88LDL	Blue				
			XL	88LDXL			USP 63.0 min		
Phlebotom	V							<u>I</u>	l
	Phlebotomy and high tactile sensitivity applications	FLEXAL™ Feel	XS	88TT20XS					VC I
			S	88TT21S	Blue	\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\			XS – L: 300 ea/bx;
			M	88TT22M					10 bx/cs
3 3 mil			171	001122101	Dide		1	I	VI.
3.3 mil			L	88TT23L) (XL: 250 ea/bx;

For more information or to order, contact your Cardinal Health sales representative.

Call 877.CARDINAL or visit cardinalhealth.com/examgloves



¹ OSHA Standard 1910.138 - Hand Protection https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.138
2 Meets the ASTM D6978 testing standard for resistance to permeation of various chemotherapy drugs. This Cardinal Health product helps healthcare facilities comply with the USP x800> standard for the safe handling of hazardous drugs.
3 This glove has been tested for puncture resistance per ASTM F1342, Standard Test Method for Protective Clothing Material Resistance to Puncture.
4 Caution: This Product (Latex Decontamination) Contains Natural Rubber Latex Which May Cause Allergic Reactions.